

SALT ANALYSIS

Analysis of anion

S.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	Colour of the salt	Blue or bluish green Green colour Light green Dark brown Pink violet Light pink or flesh colour Colourless/ white	May be copper or nickel salt Nickel salt Presence of ferrous salt (Fe^{2+}) Presence of ferric salt (Fe^{3+}) May be cobalt salt May be manganese salt Absence of Cu^{2+} , Ni^{2+} , Fe^{3+} , Mn^{2+} and cobalt
2	odour	Ammoniacal smell Rotten egg smell No characteristic smell Vinegar like smell	Presence of ammonium (NH_4^+) Presence of sulphide (S^{2-}) Absence of NH_4^+ & S^{2-}
3	Appearance of salt	crystalline powdery	May be CH_3COO^- May be sulphate, nitrate or chloride May be carbonate or sulphide
4	Solubility: a little of the substance is shaken with water	Soluble insoluble	May be sulphate, nitrate or chloride May be carbonate or sulphide
5	Flame test: salt is made to a paste with 2 drops of con. HCl in a watch glass. The paste is introduced in a non luminous flame.	Bluish green flame Pale green flame Brick red flame Crimson No characteristic coloured flame	Presence of Cu^{2+} Presence of Ba^{2+} Presence of Ca^{2+} Presence of Sr^{2+} Absence of Cu^{2+} , Ba^{2+} , Ca^{2+} , and Sr^{2+}

SYSTEMATIC ANALYSIS OF ANION

s.No.	EXPERIMENT	OBSERVATION	INFERENCE
6	Salt + dil H_2SO_4	<ol style="list-style-type: none"> 1. colourless, odourless gas (CO_2) with brisk effervescence turning lime water milky 2. colourless rotten egg smelling gas (H_2S) turning lead acetate paper black. 3. colourless gas (SO_2) turns acidified potassium dichromate paper green 4. Reddish brown gas (NO_2) turns ferrous sulphate solution black 5. no characteristic change 	<p>Presence of carbonate CO_3^{2-}</p> <p>Presence of sulphide S^{2-}</p> <p>May be SO_3^{2-} (sulphite)</p> <p>May be NO_2^-</p> <p>Absence of CO_3^{2-} and S^{2-} SO_3^{2-} and NO_2^-</p>
7.	Salt + con H_2SO_4	<ol style="list-style-type: none"> 1. colourless gas (HCl) giving white fumes with ammonium hydroxide taken at the end of the glass rod. 2. Reddish brown gas. (Br_2/NO_2) 3. deep violet vapours (I_2) 4. Colourless gas with vinegar like 	<p>Presence of chloride</p> <p>Presence of NO_3^- / Br^-</p> <p>Presence of iodide</p> <p>Presence of acetate ion</p>

		smell turns blue litmus red 5. colourless gas turns pink colour KMnO_4 into colourless 6. No characteristic change	Presence of oxalate ion Absence of Cl^- , Br^- , I^- , CH_3COO^- and oxalate ion
8.	Copper turning test: Salt + con H_2SO_4 + cu turning + heat	Reddish brown gas	Presence of Nitrate
9.	BaCl_2 test. Salt solution + barium chloride solution	No brown gas White precipitate (undissolved in dil. HCl) White PPte dissolved in con. HCl No PPte	Absence of NO_3^- Presence of SO_4^{2-} Presence of SO_3^{2-} Absence of SbO_4^{2-} & SO_3^{2-}
10.	Ammonium molybdate test: Salt + Con. HNO_3 - heated and cooled, then added to excess of Ammonium molybdate	Yellow precipitate soluble in NH_3 No yellow PPte	PO_4^{3-} is present Absence of PO_4^{3-}
For cation			
12	Nessler's reagent test Salt solution + Nessler's reagent	Reddish brown colour	NH_4^+ ion is present

Preparation of a sodium carbonate extract:

1g of the salt, 3g of sodium carbonate and 20 ml of water are taken in a beaker. Heated to boil and then filtered. The filtrate is sodium carbonate extract.

CONFIRMATORY TEST FOR ANIONS:

Confirmatory test for carbonates:

1	Salt solution + dil.HCl Pass this gas into lime water	Brisk effervescence Milky white precipitate	Carbonate is confirmed
2	Salt solution + magnesium sulphate	White precipitate	Carbonate is confirmed

Confirmatory test for sulphides

s.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	Salt + dil H_2SO_4 and lead acetate paper is dipped + dil H_2SO_4	Lead acetate paper turns black	Sulphide is confirmed
2.	Salt solution + $CdCO_3$	Yellow precipitate	Sulphide is confirmed

X Confirmatory test for sulphite:

s.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	$BaCl_2$ test. Salt solution + barium chloride solution	White precipitate (dissolved in dil. HCl)	Sulphite is confirmed
2.	Salt solution + $K_2Cr_2O_7$	Green colouration	Sulphite is confirmed

X Confirmatory test for NO_2^- :

s.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	Salt solution + dil acetic acid + $FeSO_4$ solution	A dark brown or black colouration	Nitrite is confirmed
2.	Salt solution + few drops of diPhenyamine	A deep blue colouration	NO_2^- is confirmed

Confirmatory test for -chloride:

s.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	Salt solution +dil HNO ₃ boil and cool then add AgNO ₃	Curdy white precipitate soluble in ammonium hydroxide	chloride is confirmed
2	Salt +solid MnO ₂ + con. H ₂ SO ₄	Greenish yellow gas having pungent smell.	chloride is confirmed
3	<u>Chromyl chloride test:</u> Salt + powdered Potassium dichromate + Con H ₂ SO ₄ Heat the mixture and pass red vapours evolved in to NaOH solution To the yellow solution add dil acetic acid and lead acetate.	Yellow ppte	chloride is confirmed

Confirmatory test for bromide:

s.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	Salt solution +dil HNO ₃ boil and cool then add AgNO ₃	Pale yellow precipitate	bromide is confirmed
2	Salt +solid MnO ₂ + con. H ₂ SO ₄	Reddish brown gas	bromide is confirmed

Confirmatory test for iodide:

s.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	Salt solution +dil HNO ₃ boil and cool then add AgNO ₃	yellow precipitate soluble in ammonium hydroxide	iodide is confirmed
2	Salt +solid MnO ₂ + con.	violet vapours	iodide is

	H ₂ SO ₄		confirmed
--	--------------------------------	--	-----------

Confirmatory test for nitrate:

S.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	<u>Brown ring test</u> Salt solution is treated with a freshly prepared ferrous sulphate and con; sulphuric acid is added in drops along the sides of the test tube.	Brown ring is formed	nitrate is confirmed
2	<u>Copper turning test</u> Salt is heated with cu turnings and con. Sulphuric acid.	Dark brown gas fumes of nitrogen dioxide is evolved	nitrate is confirmed

Confirmatory test for sulphate

S.No.	EXPERIMENT	OBSERVATION	INFERENCE
1	Salt solution + Barium chloride solution	A white precipitate insoluble in acids	Sulphate is confirmed.
2	Salt solution + lead acetate	A white precipitate	Sulphate is confirmed

Confirmatory test for phosphate

S.No.	EXPERIMENT	OBSERVATION	INFERENCE
-------	------------	-------------	-----------

1	Salt + Con. HNO_3 – heated and cooled, then added to excess of Ammonium molybdate	Yellow precipitate soluble in NH_3	Phosphate is confirmed
2	Salt + Magnesia mixture	White precipitate	Phosphate is confirmed

* Confirmatory test for acetate:

1	Salt + FeCl_3 solution heat for 5 min and cool	Blood colour solution	Acetate is confirmed
2	Ester test: salt + con sulphuric acid add ethyl alcohol. Warm in the water bath.	Fruity smell	Acetate is confirmed

* Confirmatory Test for oxalate:

1	Salt solution + dil acetic acid + CaCl_2 solution	White precipitate	Oxalate is confirmed
2	Dissolve the part of white ppt in dil. H_2SO_4 and add KMnO_4 solution. Warm	Pink colour disappears	Oxalate is confirmed

* Test for Nitrite

1	Salt solution + dil acetic acid. Boil to expel CO_2 add solid KI and starch solution	Deep blue colour	NO_2^- is confirmed
2			